

README file for PFMCal.

A. Butykai, P. Domínguez-García, F. M. Mor, R. Gaál, L. Forró, S. Jeney

1. ORIGINAL PFMCal.M

1.1 Software requirements

MatLab 2011a or later, Statistics toolbox. Tested in MatLab 2016a.

1.2 Installation of the software

Unzip `CalibrationSoftware.zip` archive to the destination path. No further installation steps are needed.

1.3 Executing the software

The calibration software can be executed from MatLab (MathWorks Inc.) by running the `PFMCal.m` file. Select either ‘Add to path’ or ‘Change current directory’ options.

1.4 Using the software

Please refer to Section 7 in the original manuscript: **Computer Physics Communications, 196 (2015) 599**

1.5 Test data

Test data are located in the `Testdata` folder. `Si_Spheres` archive contains experimental data measured on silica spheres. `Si_splinters_non_spherical` folder contains experimental data measured on fragmented silica particles with random shape. `Synthetic_non_spherical` folder contains 100 generated datasets with random error for the validation of the non-spherical model. `Synthetic_spherical` folder contains 100 generated datasets with random error for the validation of the spherical model.

2. EXTENDED VERSION

2.1 Software requirements

The extended code can be executed without modification using MatLab or GNU Octave. The code has been tested in Linux and Windows operating systems.

2.2 Installation of the software

Unzip `CalibrationSoftware.zip` archive to the destination path. The extended files are inside the `AUTO` folder. For running the files an installation of MatLab 2011 (or later) or GNU Octave 4.0 (or later) is needed.

Regarding Octave, the code uses `Forge*` packages `struct` and `optim`. Issues may appear in relation with those packages in Octave 4.0 (not expected in later versions). To avoid these errors, update the packages in Octave command window by using the following commands:

```
pkg install --forge struct
pkg install --forge optim
```

*<https://octave.sourceforge.io/>

2.3 Executing the software

The extended calibration software can be executed from MatLab command window or GNU Octave by running the files `PFMCal_auto.m` and `PFMCal_histo.m` located inside folder `AUTO` (Linux-style complete path: `./CalibrationSoftware/AUTO`). Therefore, for the automatic process of calibration, enter in the MatLab/Octave command window:

```
>PFMCal_auto
```

For running the thermal noise statistics method, simply write:

```
>PFMCal_histo
```

There is no need to enter a path to the data files, because it is included in those scripts. The test data are loaded from the folder `Testdata` (Path: `./CalibrationSoftware/AUTO/Testdata`). The loading of the files have been written to be compatible with the different path separators depending on the operative system (tested in Windows and Linux). By default, the water data is selected, but it can be changed to acetone data by editing directly the main files `PFMCal_auto.m` or `PFMCal_histo.m`. Inside these two files, the input data needed for each running are defined. This input data are defined in the following sections.

2.4 Using the software

2.4.1 PFMCal_auto

The analytical solution of Langevin equations for an optically trapped spherical particle with no-slip boundaries in a Newtonian fluid with hydrodynamic effects only depends numerically on the following timescales:

$$\tau_f \equiv \frac{\rho_f a^2}{\eta}; \tau_p^* \equiv \frac{m^*}{6\pi\eta a}; \tau_k \equiv \frac{6\pi\eta a}{k} \quad (1)$$

where a is the bead radius and ρ_p its density, ρ_f is the density of the fluid, while $m^* \equiv m_p + m_f/2$ is an effective mass which modifies the mass of the particle, m_p , because of the influence of hydrodynamics, m_f being the mass of the displaced fluid. By knowing the characteristics of the beads and the fluid (a , ρ_f and ρ_p), the calibration method allows to calculate these timescales and, consequently, to obtain the viscosity of the fluid, η , and the optical force spring constant, k .[†]

Additionally to these input values, we have to obtain the value of the MSD plateau which is going to be used in the computation of the timescales. The value of the MSD plateau is often difficult to obtain as data can be quite noisy in the plateau region (see Fig. 1). Errors in the estimation of this value results in the imprecise determination of η .

In order to extract a reliable plateau value, the function `getMSDPlateau` automatically calculates the second derivative of the logarithm of the MSD:

$$f(t) \equiv \frac{d^2}{dt^2} \log \text{MSD}(t)$$

[†]For further details see: M. Grimm, T. Franosch, S. Jeney, High-resolution detection of Brownian motion for quantitative optical tweezers experiments, Phys. Rev. E 86 (2012) 021912.

which is a measure of the noise affecting the MSD. The function `getMSDPlateau` further calculates a modulation function:

$$M_n \equiv \frac{\max(f_n) - \min(f_n)}{\max(f_n) + \min(f_n)}$$

with $n = 1, 2, \dots$ the n th point of $f(t)$ and determines the value of n where M_n grows abruptly above 10. The corresponding value of $f(t)$ gives a precise value of the MSD plateau and is used, together with the input data a , ρ , ρ_p and T to automatically calculate the calibration curves for the MSD and VAF, and to obtain β_{MSD} , β_{VAF} , k and η . Fig. 1 shows $f(t)$ and M_n calculated from the acetone MSD data contained in the `AUTO/Testdata` folder.[‡]




Figure 1. Left: Logarithm of the second derivative of the mean-squared displacement ($f(t)$ in the text), showing the noise affecting the calculation of the plateau of the MSD. The red circle marks the approximate point chosen by the user where the noise begins to be relevant. Right: Modulation function of $f(t)$ (M_n in the text): the beginning of noise is the first peak in the figure.

Regarding the PSD, we use the classic methodology for Newtonian fluids, by fitting a Lorentzian curve to the measured data to obtain its corner frequency, f_c , from where the corresponding calibration factor β_{PSD} and k can be obtained. The PSD of the bead in the fluid has a Lorentzian form $\text{PSD}(f) = S_0 f_c^2 / (f_c^2 + f^2)$, where S_0 is a constant. This quantity can be approximated as $\text{PSD}(f) \simeq S_0 f_c^2 f^{-2}$ in the limit of high frequencies. Then, the calibration factor is retrieved using $\beta_{\text{PSD}} = k_B T / (6\pi^3 \eta a S_0 f_c^2)$. Note here that η and a have to be previously known to obtain this β_{PSD} , and therefore we use the η value obtained from the calibration of the MSD and VAF.

The execution of `PFMCal_auto.m` finally provides Fig. 2, which plots measured MSD, VAF and PSD together with the fitted model functions. The characteristic timescales obtained by the fitting are displayed in the graphs with vertical lines. This graph is automatically saved in a folder created by the software named `output`. In this folder we save the calibration graph in ‘pdf’ and ‘eps’ formats and a text file named `output_data.txt` which contains all the input and output data.

[‡]This calculation is not shown during the execution of the code but these graphs can be shown and saved by changing a boolean parameter in `getMSDPlateau` function.




Figure 2. Experimental MSD and the VAF with their theoretical curves. The experimental PSD is indicated by a fitting to the Lorentzian curve. Vertical lines: t_f corresponds to τ_f , t_{peff} is τ_p^* , t_k is τ_k . The time window in which the described methods can be applied are delimited by t_{fit1} and t_{fit2} .

2.4.2 PFMCAl.histo.m

Additionally to the methods explained above, we also include a calibration method by using the experimental histogram of particle positions. By calculating the probability density $p(x)$, the potential can be deduced as $E(x) = -k_B T \ln p(x) + C$, where C is a constant related to the potential offset and can be neglected. In the case of a spherical bead in an ideal optical trap, $E(x) = \frac{1}{2}k_E x^2$ in a harmonic approximation, only depending on the trap stiffness, k . Here, $p(x)$ is a Gaussian distribution and, through the equipartition theorem: $\frac{1}{2}k \langle x^2 \rangle = \frac{1}{2}k_B T$, we have $k_{\sigma^2} \equiv k_B T / \sigma^2$, where $\sigma^2 = \langle x^2 \rangle$ is the variance of the probability density.

If we use some additional method to obtain a value for the trap stiffness, i.e., by calculating averaged values from the calibration of MSDs for independent water measurements, we can calculate two additional calibration factors, β_E , from the harmonic potential curve and β_{σ^2} , from the Gaussian distributions.




Figure 3. Calibration factors obtained from the fitting to a harmonic potential curve, β_E , and from the Gaussian distribution, β_{σ^2} .

The fittings and calculation of these calibration factors is done by using the script `PFMCAl_histo.m`. The histograms data are read from the `AUTO/Testdata` folder and the k value is defined inside `PFMCAl_histo.m`. This code fits the data to a Gaussian function and to a harmonic potential, obtaining calibration values very similar to the ones determined with the calibration methodology of `PFMCAl.m` and `PFMCAl_auto.m` (see Table in the new paper submission). The fitting graph is like Fig. 3, which is saved in the `output` folder. The input and output data regarding this calibration code are saved in a file named `output_data_histo.txt`.

All the fittings results (texts and graphs for water and acetone data) are shown in the folder `./AUTO/results_auto`.